

CHAPTER 11

PORTABLE BREATHING OXYGEN REGULATOR AND CYLINDER

TYPE MA-1 (CLIFTON PRECISION TYPE)


Section 11-1. Description

11-1. GENERAL.

11-2. The Portable Breathing Oxygen Regulator and Cylinder, Type MA-1, using regulator P/N 3260007-0101, is manufactured by Litton Life Support formerly Clifton Precision (CAGE 99251) (figure 11-1). The regulator and cylinder, is a self-contained portable breathing device capable of supplying breathing oxygen to flight personnel for normal or emergency use. The regulator type A-21 (P/N 3260007-0101 and P/N 3341007-1), hereafter referred to as the A-21 which forms part of the unit, is a demand/pressure breathing type regulator which will deliver oxygen to the user upon demand or provide a positive pressure to the mask depending upon the positioning of the selector knob. The unit is intended for use with the full face smoke mask (P/N 651-280) or the MBU-12 oxygen mask configuration. During normal operation the selector knob is positioned in the norm position and deliver 100% oxygen upon demand. When the selector knob is placed in the 30M, 42M, or EMER position, the unit will deliver 100% oxygen at a positive pressure of 1.6 to 14.0 inH₂O depending upon the positioning of the selector knob from sea level to the service ceiling of the aircraft. See table 11-1 for leading particulars.

11-3. CONFIGURATION.

11-4. The regulators are supplied in one basic configuration: low pressure (50 to 500 psig operating pressure range).



011001

Figure 11-1. Portable Breathing Oxygen Cylinder and Regulator, Type MA-1

Table 11-1. Leading Particulars

CHARACTERISTICS	
A-21	P/Ns 3260007-0101 and 3341007-1 50 to 500 psig
Mounting	MA-1 Portable walk around cylinder
Operating Altitude	Range Sea Level to 43,000 ft
Visual Indicator	Pressure gage
Selector knob:	
NORM	100% Oxygen on demand
30M	Provides 100% oxygen at a positive pressure of 1.6 to 2.4 inH ₂ O
42M	Provides 100% oxygen at a positive pressure of 5.5 to 6.5 inH ₂ O
EMER	Provides 100% oxygen at a positive pressure of 12.0 to 14.0 inH ₂ O
Overall Dimensions:	
Pressure gage	
Length	5 7/8 in.
Width	5 1/4 in.
Height	3 5/8 in.
Weight	2.2 lbs.

11-5. SERVICE LIFE.

NOTE

Low Pressure Oxygen Cylinders, Type MS21227-1 used on MA-1 Portable Emergency Oxygen System, do not have a Department of Transportation (D.O.T.) ICC number

permanently stamped in the neck of the cylinder and therefore do not require hydrostatic testing. These cylinders are painted yellow in accordance with MIL-STD-101.

11-6. The regulator shall remain in service as long as it continues to function correctly and does not require other than authorized adjustment or repair.

Section 11-2. Modifications

11-7. GENERAL.

11-8. There are no authorized modifications to the Type A-21 regulator.

Section 11-3. Performance Test Sheet Preparation

11-9. GENERAL.

11-10. Preparation of the oxygen regulator Performance Test Sheet requires that, through the use of various graphs, actual flows given in applicable directives and provided in this section be converted to indicated flows.

11-11. Flows provided in applicable directives are stated in liters per minute (lpm) and are not measurable by the manometers used in oxygen regulator test stands. The flows must be converted to inches of water pressure (inH₂O), the form of measurement which can be read on the test stand manometers.

11-2 Change 4

NOTE

The various graphs supplied with each oxygen system components test stand, models 1172AS100 and 1316AS100 are used in converting flows. The graphs supplied are not interchangeable between test stands.

11-12. The information provided in the tables in this section are to be recorded on the Performance Test Sheet (figure 11-2).

11-13. The Performance Test Sheet (figure 11-2) is a sample only, but may be reproduced for local use.

11-14. The following tests require conversion of flows from actual lpm to indicated inH₂O.

1. Flow Suction Test.
2. Pressure Breathing Test.

11-15. REGULATOR PERFORMANCE TESTS.

11-16. FLOW SUCTION TEST. The flow suction test shall be performed at sea level with the selector knob in the norm position. Actual (lpm) flows are converted to indicated (inH₂O) output flows by using the sea level output graph N₂ line. Convert the actual flows as follows:

NOTE

Test stand output flow graphs may vary in makeup, according to the activity performing the test stand calibration. Some test stands may have a single output graph with various altitude lines, while others may have separate graphs for each altitude. Ensure specified graph is used.

1. Locate the desired lpm line (figure 11-2) at the bottom of the sea level output graph.

2. Trace selected lpm line up to where it intersects the N₂ line.

3. Trace the line from point of intersection across the graph to the left-hand column to determine indicated output inH₂O.

4. Enter this figure in the appropriate block on the Performance Test Sheet.

5. Repeat steps 1 through 4 for all output flows (LPM) given in figure 11-2.

11-17. PRESSURE BREATHING TEST. The Pressure Breathing Test shall be performed at sea level with the selector knob in the 30M, 42M, and EMER positions. Actual (lpm) flows are converted to indicated (inH₂O) output flows by using the sea level output graph N₂ line. Convert the actual flows as follows:

NOTE

Test stand output flow graphs may vary in makeup, according to the activity performing the test calibration. Some test stands may have a single output graph with various altitude lines, while others may have separate graphs for each altitude. Ensure specified graph is used.

1. Locate the desired lpm line (figure 11-2) at the bottom of the sea level output graph.

2. Trace selected lpm line up to where it intersects the N₂ line.

3. Trace the line from point of intersection across the graph to the left-hand column to determine indicated output inH₂O.

4. Enter this figure in the appropriate block on the Performance Test Sheet.

5. Repeat steps 1 through 4 for all output flows (lpm) given in figure 11-2.

Section 11-4. Maintenance**11-18. GENERAL.**

11-19. Maintenance of the Portable Breathing Oxygen Cylinder and Regulator, Type MA-1, is limited to

inspection, testing, adjustment of the pressure breathing assembly, and tightening of loose screws and fittings. Ensure that anti-seize tape is used on pipe threads. This section contains only the procedural steps necessary to meet these requirements.

NAVAIR 13-1-6.4-2

REGULATOR PERFORMANCE TEST SHEET
TYPE A-21
3260007-0101 and 3341007-1

DATE: _____ TYPE: _____ SERIAL NUMBER: _____

TEST STAND SERIAL NUMBER: _____ TESTED BY: _____ CDI: _____

1. OVERALL LEAKAGE TEST: _____

2. PRESSURE GAGE SCALE AND ERROR TEST:

INLET PRESSURE	TOLERANCE	PRESSURE GAGE READING
200 PSIG	± 25 PSIG	
300 PSIG	± 25 PSIG	
400 PSIG	± 25 PSIG	

3. FLOW SUCTION TEST:

SELECTOR KNOB SETTING	INLET PRESSURE	ACTUAL OUTPUT LPM	INDICATED OUTPUT INH ₂ O	READING	MAXIMUM ALLOWABLE SUCTION
NORM	50	10			-0.4
NORM	50	50			-1.4
NORM	150	100			-1.4
NORM	500	10			-1.4
NORM	500	50			-1.4
NORM	500	100			-1.4

Figure 11-2. Performance Test Sheet Preparation (Sheet 1 of 2)

4. PRESSURE BREATHING TEST:

SELECTOR KNOB SETTING	INLET PRESSURE (PSIG)	ACTUAL OUTPUT LPM	INDICATED OUTPUT IN H ₂ O	READING	MAXIMUM ALLOWABLE SUCTION
30M	150	10			+1.6 TO +2.4
30M	150	90			-1.0 NOTE 1
30M	150	0	0		+1.0 NTOE 2
42M	150	10			+5.5 TO +6.5
42M	150	90			-1.5 NOTE 1
42M	150	0	0		+1.5 NOTE 2
EMER	150	10			+12.0 TO +14.0
EMER	150	90			-2.0 NOTE 1
EMER	150	0	0		+2.0 NOTE 2

5. OXYGEN CYLINDER AND REGULATOR PURGE _____

FILL PRESS.	DEplete	FILL PRESS.	DEplete	FILL PRESS.	DEplete
450 PSIG		450 PSIG		450 PSIG	

NOTES:

1. MAXIMUM ALLOWABLE PRESSURE DECREASE FROM 10 LPM READING.
2. MAXIMUM ALLOWABLE PRESSURE INCREASE ADDED TO 10 LPM READING.

Figure 11-2. Performance Test Sheet Preparation (Sheet 2 of 2)

NOTE

Upon completion of any maintenance action (e.g., inspection, adjustment, repair, etc), be sure to complete the required maintenance data collection system form.

11-20. INSPECTIONS.

11-21. If the regulator fails any inspection or testing requirements and cannot be adjusted within limits which are authorized, it shall be removed from service and sent to Air Force Depot Level Maintenance for repair and replaced with a Ready-For-Issue (RFI) component.

11-22. Procedural steps outlined in this section are listed under the inspection cycle in which they are required, and in the sequence in which they normally occur.

11-23. TURNAROUND/PREFLIGHT/POST-FLIGHT/TRANSFER INSPECTIONS. The Turn-around/Preflight/Postflight/Transfer Inspections consist of a visual inspection performed in conjunction with the aircraft inspection requirements for the aircraft in which the regulators are installed. See [table 11-2](#) for assistance in troubleshooting. To perform the inspection, visually inspect the following:

1. Legibility of all markings.
2. Low or improper reading on regulator pressure gage.
3. Selector knob in NORM position.
4. Regulator and surrounding area for freedom from dirt and hydrocarbons.
5. Delivery hose and connector for cuts, fraying, kinking, hydrocarbons and general condition.

11-24. If discrepancies are found or suspected, Maintenance Control shall be notified.

11-25. Regulators which do not pass inspection and cannot be repaired in the aircraft, shall be removed and replaced by Ready-For-Issue (RFI) regulators. Non-RFI regulators shall be forwarded to the nearest maintenance activity having repair capability.

11-26. ACCEPTANCE/SPECIAL/DAILY INSPECTIONS. The Acceptance/Special/Daily Inspections consist of a visual inspection followed by a Functional Test. These inspections and tests shall be performed in

conjunction with the aircraft inspection requirements for the aircraft in which the regulators are installed. Refer to [table 11-2](#) for assistance in troubleshooting. To perform the inspection, proceed as follows:



When working with oxygen, make certain that clothing, tubing fittings, and equipment are free of oil, grease, fuel, hydraulic fluid, or any combustible liquid. Fire or explosion may result when even slight traces of combustible material come in contact with oxygen under pressure.

11-27. Visually inspect the regulators in accordance with [paragraph 11-23](#).

11-28. FUNCTIONAL TEST. To perform the Functional Test, proceed as follows:

1. Place selector knob in NORM position.
2. Connect oxygen hose to quick-disconnect, place mask to face and inhale, there should be no restriction inhaling or exhaling.
3. Turn selector knob to 30M, 42M, and EMER positions. There should be a progressively higher flow rate when moving selector knob from 30M to 42M and from 42M to EMER positions.

11-29. Upon completion of Functional Test, secure regulator as follows:

1. Place selector knob in NORM position.
2. Disconnect mask hose from regulator.

11-30. If discrepancies are found or suspected, maintenance control shall be notified.

11-31. Regulators which do not pass inspection and cannot be repaired in the aircraft, shall be removed and replaced by RFI regulators. Non-RFI regulators shall be forwarded to the nearest maintenance activity having repair capability.

11-32. CALENDAR/PHASED/SDLM INSPECTIONS. Calendar, Phased or SDLM Inspections require removal of the regulators from the aircraft. See applicable planned maintenance system (PMS) publications for specified intervals. In no case shall the interval exceed 448 days. Upon removal from the aircraft, regulators shall be visually inspected ([paragraph 11-23](#)) and bench tested ([paragraph 11-34](#)).

Table 11-2. Troubleshooting (Daily, Preflight, Special, Turnaround, Transfer, and Acceptance Inspections)

Trouble	Probable Cause	Remedy
Oxygen cylinder pressure gage fails to indicate proper pressure.	Defective gage.	Replace regulator.
	Low cylinder pressure.	Refill.
	Defective inlet assembly.	Replace regulator.
Oxygen not available at mask with selector knob in NORM position with proper pressure source to regulator.	Hose to mask is kinked.	Straighten hose and reposition outlet.
	Regulator not functioning properly.	Replace regulator.
Oxygen not available at mask with proper source to regulator with selector knob in 30M, 42M, and EMER positions.	Kinked or other malfunction between hose and mask.	Replace or readjust equipment as necessary.
	Regulator malfunctioning or out of adjustment.	Replace regulator.
Gage pressure drops when regulator is not in use.	Loose or leaking connections.	Tighten or replace connections as necessary.
	Defective regulator.	Replace regulator.

11-33. Aircraft type A-21 regulators failing the Bench Test and cannot be adjusted within limits shall be shipped to Depot Level repair.

changes. Troubleshooting tables are provided following each test.

11-34. BENCH TEST.



Because of possible vacuum pump explosion, only water-pumped nitrogen, Type I, Class I, Grade B (Fed Spec BB-N-411) shall be used in testing oxygen regulators.

11-35. Bench Test shall be performed on aircraft oxygen regulators prior to being placed in service, and during the Phase/Calendar or SDLM (Standard Depot Level Maintenance) Inspection cycle of the aircraft in which installed. See applicable planned maintenance system (PMS) publications for specific intervals. In no case shall the interval exceed 448 days. The regulators shall also be subjected to a Bench Test if malfunction is suspected, and after repair or replacement of damaged parts.

For oxygen test stands and purging equipment, use only nitrogen from gray cylinders marked NITROGEN OIL FREE in white letters. Two 3-inch wide black bands mark the tops of these cylinders. Do not use 3500 psig nitrogen cylinders. These cylinders cannot be certified contaminant free.

11-36. Bench Test shall be performed using Oxygen System Components Test Stand, Model 1172AS100 or 1316AS100. Refer to the appropriate ground support equipment manual for identification of test stand controls and indicators referred to in Bench Test.

NOTE

Tests are arranged so they proceed from one test to the next with a minimum of flow

11-37. Due to the complexity of the 1172AS100 and 1316AS100 test stand, it is essential that the operator become thoroughly familiar with the test stand prior to performing Bench Test. Refer to appropriate ground support equipment manual.

11-38. DEPLETION OF OXYGEN CYLINDER AND REMOVAL OF OXYGEN REGULATOR FROM CYLINDER. To deplete oxygen and remove regulator from cylinder, proceed as follows:

WARNING

When working with oxygen, make certain that clothing, tubing fittings, and equipment are free of oil, grease, fuel, hydraulic fluid, or any combustible liquid. Fire or explosion may result when even slight traces of combustible material come in contact with oxygen under pressure.

Deplete oxygen from the oxygen cylinder in an area free from dirt, grease, oil, or other hydrocarbons.

1. Turn selector knob to the EMER position and completely deplete oxygen from oxygen cylinder.

CAUTION

Use care not to damage regulator or cylinder during regulator removal from cylinder. Make certain wrench is firmly on hexagon flats at base of regulator.

NOTE

For step 2 index numbers for removing regulator from oxygen cylinder, refer to [figure 11-3](#).

2. Remove regulator (1) from cylinder (4) by firmly grasping adapters (2) at base of regulator and hexagon boss (3) at cylinder neck with two 3/4 inch wrenches.

3. Cap open port of cylinder with plug and stow cylinder in bin free from dirt, grease, oil, and other hydrocarbons.

11-39. MODIFICATION OF OXYGEN HOSE COUPLING ASSEMBLY. To modify the oxygen hose coupling assembly (MS22016), proceed as follows:

1. Modify oxygen hose coupling assembly by removing spring and valve assembly and cross bar support from inside of oxygen hose coupling assembly; to provide a clear passage with no restrictions for flows being drawn.

2. Attach modified oxygen hose coupling to piezometer (26) using a short hose maintaining one inch between the oxygen hose coupling and piezometer (26).

WARNING

Ensure altitude chamber is configured in accordance with NAVAIR 17-15BC-21, WP003 00, Figure 3 sheets 2 thru 4 as applicable. Ensure High Pressure or Low Pressure Hose Assembly listed in NAVAIR 17-15BC-21, WP031 00, Figure 1 or Figure 2 is attached to N₂ Input Connection (18) or Tee Connection (28) in altitude chamber as applicable for the oxygen regulator being tested. Remove hose assembly not being used and cap connection (18) or (28) when not in use. For regulators requiring inlet pressures greater than 175 psig, the High Pressure Hose Assembly in NAVAIR 17-15BC-21, WP031 00, Figure 1 shall be used.

11-40. OVERALL LEAKAGE TEST. To perform the Overall Leakage Test, proceed as follows:

Materials Required		
Quantity	Description	Reference Number
As Required	Nitrogen, Oil-free, Water Pumped, Type I, Class I, Grade B	Fed Spec BB-N-411 NIIN 00-985-7275
As Required	Tape, Anti-seize	MIL-T-27730A, NIIN 00-889-3535 (CAGE 81349)
As Required	Compound, Leak Detection, Type 1	MIL-L-25567

Support Equipment Required		
Quantity	Description	Reference Number
1	Adapter, Straight	MS24398D8-4 (CAGE 96906) NIIN 00-684-7292
1	Oxygen Systems Components Test Stand	1172AS100 or 1316AS100

1. Apply anti-seize tape to inlet fitting of the regulator.
2. Attach adapter to inlet fitting of the regulator.
3. Ensure all test stand valves are secure and open N₂ supply cylinder.
4. Attach inlet of regulator to N₂ INPUT connection (18) in the altitude chamber.

5. Turn INLET PRESSURE ON/OFF valve (L) ON.
6. Place selector knob of regulator in NORM position.
7. Using high pressure regulator (Q) slowly apply 500 psig to regulator.
8. Using leak detection compound, draw a film over outlet of regulator. No leakage allowed.
9. Using leak detection compound, check for leakage on pressure gage, relief valve, inlet fitting, and filler valve. No leakage allowed.
10. Turn INLET PRESS. ON/OFF valve (L) OFF.
11. If leakage is excessive, locate probable cause using troubleshooting chart (table 11-3).
12. Turn high pressure regulator (Q) to vent.
13. Bleed pressure from regulator by placing selector knob in 30M position. Bleed pressure from test stand with system bleed valve (S).
14. Turn system bleed valve (S) OFF and place regulator selector knob in NORM position.

11-41. PRESSURE GAGE SCALE AND ERROR TEST. To perform the Pressure Gage Scale and Error Test, proceed as follows:

Materials Required		
Quantity	Description	Reference Number
As Required	Nitrogen, Oil-free, Water Pumped, Type I, Class I, Grade B	Fed Spec BB-N-411 NIIN 00-985-7275

Support Equipment Required		
Quantity	Description	Reference Number
1	Oxygen Systems Components Test Stand	1172AS100 or 1316AS100

1. Turn INLET PRESSURE ON/OFF valve (L) ON.



When applying inlet pressure using high pressure regulator (Q), apply pressures slowly.

2. Using high pressure regulator (Q), slowly increase pressure to each test pressure specified on Performance Test Sheet (figure 11-2).
3. Check tolerance by comparing regulator pressure gage reading with test stand regulated high pressure gage (10).
4. Regulators failing the pressure gage scale and error test shall be BCM to Air Force Depot Level Maintenance for repair.
5. Turn high pressure regulator (Q) to vent position.
6. Bleed test stand using system bleed valve (S). Bleed pressure from regulator by placing selector knob in 30M position. Close system bleed valve (S).

11-42. FLOW SUCTION TEST. To perform the Flow Suction Test, proceed as follows:

Materials Required		
Quantity	Description	Reference Number
As Required	Nitrogen, Oil-free, Water Pumped, Type I, Class I, Grade B	Fed Spec BB-N-411 NIIN 00-985-7275

Support Equipment Required		
Quantity	Description	Reference Number
1	Oxygen Systems Components Test Stand	1172AS100 or 1316AS100




Figure 11-3. Type MA-1 Portable Breathing Oxygen Cylinder and Regulator

011003



LOW PRESSURE REGULATOR (N) can only be used when applying pressures below gage guard setting (165 to 175 psig) to an item under test. For pressure above gage guard setting, HIGH PRESSURE REGULATOR (Q) must be used.

1. Place selector knob in NORM position.
2. Attach regulator outlet to piezometer (26).
3. Turn on vacuum pump.
4. Using LOW PRESS. REGULATOR (N), set inlet pressures of 50 and 150 psig as specified on Performance Test Sheet (figure 11-2).

NOTE

When drawing flows with OUTPUT valve (C), ensure inlet pressure is maintained.

5. Place PRESSURE SELECTOR valve (D) in H₂O position.
6. Using OUTPUT valve (C), set flows specified on performance test sheet on output manometer (1). Suction values will be displayed on PRESS./SUC-TION manometer (4). Record suction readings on Performance Test Sheet.
7. Close OUTPUT valve (C) and back out on low press. regulator (N).
8. Using HIGH PRESSURE REGULATOR (Q), slowly set inlet pressure of 500 psig as specified on Performance Test Sheet.
9. Using OUTPUT valve (C), set flows specified in Performance Test Sheet on output manometer (1). Suction values will be displayed on press/suction manometer (4). Record suction readings on Performance Test Sheet.

10. Regulators failing flow suction test shall be BCM to Air Force Depot Level Maintenance for repair.

11. Close OUTPUT valve (C).
12. Turn HIGH PRESSURE REGULATOR (Q) to vent.
13. Bleed pressure from test stand by opening system bleed valve (S). Bleed pressure from regulator by turning selector knob to 30M position. Close system bleed valve (S).

11-43. PRESSURE BREATHING TEST. To perform Pressure Breathing Test, proceed as follows:

Materials Required

Quantity	Description	Reference Number
As Required	Nitrogen, Oil-free, Water Pumped, Type I, Class I, Grade B	Fed Spec BB-N-411 NIIN 00-985-7275

Support Equipment Required

Quantity	Description	Reference Number
1	Oxygen Systems Components Test Stand	1172AS100 or 1316AS100

1. Place selector knob in 30M position.
2. Using LOW PRESS. REGULATOR (N), set inlet pressure of 150 psig as specified on Performance Test Sheet.

NOTE

When drawing flows with output valve (C), ensure inlet pressure is maintained.

3. Using OUTPUT valve (C), draw a 10 lpm flow as indicated on output manometer (1).

Table 11-3. Troubleshooting (Overall Leakage Test)

Trouble	Probable Cause	Remedy
Inlet fitting leaking.	Loose inlet fitting.	Tighten inlet fitting.
	Damaged threads on inlet fitting.	BCM regulator.
Pressure gage leaking.	Loose pressure gage.	Tighten pressure gage.
	Damage pressure gage.	BCM regulator.
Oxygen filler valve leaking.	Defective filler valve.	BCM regulator.
Relief valve leaking.	Defective relief valve.	BCM regulator.
Leaking first stage bellows assembly.	Defective first stage bellows or gasket.	BCM regulator.
Leakage at demand valve port.	Damaged tilt valve.	BCM regulator.

4. Reading on PRESS./SUCTION manometer (4) shall be within minimum and maximum allowable inH₂O specified on Performance Test Sheet. Record reading on Performance Test Sheet.

NOTE

The initial value obtained at a 10 lpm flow at 30M, 42M, and EMER selector knob settings establishes the base line for determining minimum allowable reading for 90 lpm flows and maximum allowable reading for 0 lpm flow.

Example: Computing minimum allowable reading for 90 lpm at 50 psig, 30M selector knob setting.

- A. 10 lpm reading at 50 psig = 1.8 inH₂O
- B. Minus maximum allowable pressure drop = 1.0 inH₂O
- C. Minimum allowable reading for 90 lpm flow at 50 psig = .8 inH₂O

Example: Computing maximum allowable reading for 0 lpm flow at 50 psig, 30M selector knob setting.

- A. 10 lpm reading at 50 psig = 1.8 inH₂O

- B. Plus maximum allowable pressure increase = 1.0 inH₂O

- C. Maximum allowable reading for 0 lpm flow at 50 psig = 2.8 inH₂O

5. Using output valve (C), draw a 90 lpm flow as indicated on output manometer (1).

6. Reading on PRESS./SUCTION manometer (4) shall not decrease more than the specified inH₂O allowable decrease on the Performance Test Sheet. Record reading on Performance Test Sheet.

7. Close output valve (C).

8. Reading on PRESS./SUCTION manometer (4) shall not increase more than the specified inH₂O allowable increase on the Performance Test Sheet. Record reading on Performance Test Sheet.

9. Repeat [steps 2 through 8](#) with selector knob in the 42M and EMER positions.

10. If the regulator fails the Pressure Breathing Test, refer to [paragraph 11-44](#).

11. Turn off vacuum pump. Close N₂ supply cylinder. Using LOW PRESS. REGULATOR (N) and system bleed valves (S), relieve all pressure in test stand. Secure all test stand valves and disconnect A-21 regulator from test stand.

11-44. ADJUSTMENT PROCEDURES FOR 30M, 42M, AND EMER. SPRING GUIDE SCREWS. To adjust 30M, 42M, and EMER. Spring Guide Screws, proceed as follows:

Materials Required

Quantity	Description	Reference Number
As Required	Nitrogen, Oil-free, Water Pumped, Type I, Class I, Grade B	Fed Spec BB-N-411 NIIN 00-985-7275

Support Equipment Required

Quantity	Description	Reference Number
1	Oxygen Systems Components Test Stand	1172AS100 or 1316AS100

NOTE

To adjust spring guide screws, refer to [figure 11-4](#) for steps 1 through 17.

1. Remove cover (2) from regulator by removing four screws (1).
2. Ensure all test stand valves are secure, then open N₂ supply cylinder.
3. Connect inlet of regulator to N₂ INPUT connection (18).
4. Connect outlet of regulator to piezometer (26).
5. Turn on vacuum pump.
6. Turn INLET PRESSURE ON/OFF valve (L) ON.
7. Using LOW PRESS. REGULATOR (N), apply 150 psig to regulator as indicated on regulated low press. gage (11).
8. Open OUTPUT valve (C) and draw a 10 lpm flow as indicated on output manometer (1).
9. Turn PRESS. selector valve (D) to H₂O position.
10. Turn regulator selector knob clockwise to 30M position.

NOTE

When performing [steps 12, 14, and 16](#) to increase outlet pressure, turn spring guide screws (3, 4, and 5) clockwise. To decrease outlet pressure, turn spring guide screws (3, 4, and 5) counterclockwise.

11. With 10 lpm flow indicated on output manometer (1), adjust spring guide screw (3) to obtain a reading of 1.6 to 2.4 inH₂O as indicated on PRESS./SUCTION manometer (4).

12. Turn regulator selector knob clockwise to 42M position.

13. With 10 lpm flow indicated on output manometer (1), adjust spring guide screw (4) to obtain a reading of 5.5 to 6.5 inH₂O as indicated on PRESS./SUCTION manometer (4).

14. Turn regulator selector knob clockwise to EMER position.

15. With 10 lpm flow indicated on output manometer (1), adjust spring guide screw (5) to obtain a reading of 12.0 to 14.0 inH₂O as indicated on PRESS./SUCTION manometer (4). Close OUTPUT valve (C). Place selector knob in NORM position.

16. Install cover (2) on regulator and secure in place with four screws (1).

NOTE

At this time, it is necessary to perform the Pressure Breathing Test ([paragraph 11-43](#)) in its entirety.

11-45. OXYGEN CYLINDER AND REGULATOR PURGE. To perform purge, proceed as follows:

Materials Required

Quantity	Description	Reference Number
As Required	Aviator's Breathing Oxygen, Type I	MIL-O-27210
As Required	Tape, Anti-seize	MIL-T-27730A, NIIN 00-889-3535 (CAGE 81349)

NAVAIR 13-1-6.4-2

WARNING

Do not use oxygen system components test stand to regulate the oxygen purge pressure.

1. Using anti-seize tape, tape the pipe threads of the inlet adapter.

WARNING

Use extreme care not to damage cylinder and regulator threads when attaching regulator to oxygen cylinder.

2. Install oxygen regulator on oxygen cylinder.
3. Connect regulator filler adapter to a regulated source of aviator's breathing oxygen.
4. Slowly charge the oxygen cylinder to 450 psig.

5. Shut off oxygen source and disconnect regulator and oxygen cylinder from regulated source of aviator's breathing oxygen.

6. Turn selector knob to 30M position and oxygen from oxygen cylinder through regulator.

7. Turn selector knob to NORM position and repeat [steps 3 through 6](#), twice.

NOTE

After purging is completed, fill MA-1 portable walk around unit to 450 psig.

All equipment forwarded from the Organizational Level Maintenance to the AIMD shall be accompanied by the appropriate forms in accordance with OPNAVINST 4790.2 Series. The test stand operator and CDI shall sign the Performance Test Sheet, and the original or a copy shall be forwarded to the organizational custodian. Upon completion of the Bench Test or Calendar Inspection, the organizational custodian shall retain the appropriate forms in accordance with OPNAVINST 4790.2 Series.




Figure 11-4. Adjustment of 30M, 42M and EMER Springs

011004

THIS PAGE INTENTIONALLY LEFT BLANK.